

ARG56999 anti-Influenza A Virus H5N1 Hemagglutinin antibody [2B7]

Package: 50 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [2B7] recognizes Influenza A Virus H5N1 Hemagglutinin
Tested Reactivity	Influenza A
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Clone	2B7
Isotype	IgG1, kappa
Target Name	Influenza A Virus H5N1 Hemagglutinin
Species	Influenza A virus
Immunogen	Recombinant fragment around aa. 17-338 of Influenza A Virus H5N1 Hemagglutinin.
Conjugation	Un-conjugated


Application Instructions

Application table	Application	Dilution
	WB	1:3000

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.4), 0.02% Sodium azide and 10% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	10% Glycerol
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG56999 anti-Influenza A Virus H5N1 Hemagglutinin antibody [2B7] WB image

Western blot: 50 ng of 1) H1N1, 2) H3N2 and, 3) H5N1 HA1 recombinant protein stained with ARG56999 anti-Influenza A Virus H5N1 Hemagglutinin antibody [2B7] at 1:3000.
