

**ARG57343**  
**anti-BCAP29 antibody**Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes BCAP29
Tested Reactivity	Hu, Ms, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	BCAP29
Species	Human
Immunogen	Recombinant Protein of Human BCAP29.
Conjugation	Un-conjugated
Alternate Names	BAP29; B-cell receptor-associated protein 29; Bap29; BCR-associated protein 29

### Application Instructions

Application table	Application	Dilution
	IHC-P	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	SW480	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.


## Bioinformation

---

Gene Symbol	BCAP29
Gene Full Name	B-cell receptor-associated protein 29
Function	May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi. May be involved in CASP8-mediated apoptosis (By similarity). [UniProt]
Calculated Mw	28 kDa


## Images

---



ARG57343 anti-BCAP29 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Rat kidney tissue stained with ARG57343 anti-BCAP29 antibody at 1:100 dilution.



SW480

ARG57343 anti-BCAP29 antibody WB image

Western blot: SW480 cell lysate stained with ARG57343 anti-BCAP29 antibody.