

## ARG59146 anti-Testosterone antibody [4E1G2]

Package: 250 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [4E1G2] recognizes Testosterone
Tested Reactivity	Other
Tested Application	IHC-P, RIA
Specificity	In competitive binding immunoassay, this antibody reacts with testosterone 100%, 11-beta-hydroxy testosterone 3.3%.
Host	Mouse
Clonality	Monoclonal
Clone	4E1G2
Isotype	IgG1, kappa
Target Name	Testosterone
Species	Others
Immunogen	BSA-conjugated Testosterone-3-CMO.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	IHC-P	0.25 - 0.5 µg/ml
	RIA	Assay-dependent
Application Note	IHC-P: Antigen Retrieval: Boil tissue section in 10 mM Citrate buffer (pH 6.0) for 10-20 min, followed by cooling at RT for 20 min. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS, 0.05% Sodium azide and 0.1 mg/ml BSA.
Preservative	0.05% Sodium azide
Stabilizer	0.1 mg/ml BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed


before use.

**Note**

For laboratory research only, not for drug, diagnostic or other use.

**Images**

---



**ARG59146 anti-Testosterone antibody [4E1G2] IHC-P image**

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human testicular carcinoma stained with ARG59146 anti-Testosterone antibody [4E1G2]. Antigen Retrieval: Boil tissue section in 10mM Citrate buffer (pH 6.0) for 10-20 min followed by cooling at RT for 20 min.