

ARG59374 anti-THRSP antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes THRSP
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	THRSP
Species	Human
Immunogen	KLH-conjugated synthetic peptide between aa. 18-48 of Human THRSP.
Conjugation	Un-conjugated
Alternate Names	Lpgp; SPOT14; Thyroid hormone-inducible hepatic protein; S14; LPGP1; Spot 14 protein; THRP

Application Instructions

Application table	Application	Dilution
	WB	1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	SK-BR-3	
Observed Size	~ 18 kDa	


Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide.
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	THRSP
Gene Full Name	thyroid hormone responsive
Background	The protein encoded by this gene is similar to the gene product of S14, a rat gene whose expression is limited to liver and adipose tissue and is controlled by nutritional and hormonal factors. This gene has been shown to be expressed in liver and adipocytes, particularly in lipomatous modules. It is also found to be expressed in lipogenic breast cancers, which suggests a role in controlling tumor lipid metabolism. [provided by RefSeq, Jul 2008]
Function	Plays a role in the regulation of lipogenesis, especially in lactating mammary gland. Important for the biosynthesis of triglycerides with medium-length fatty acid chains. May modulate lipogenesis by interacting with MID1IP1 and preventing its interaction with ACACA (By similarity). May function as transcriptional coactivator. May modulate the transcription factor activity of THRβ. [UniProt]
Calculated Mw	17 kDa
Cellular Localization	Nucleus. Cytoplasm. [UniProt]

Images



ARG59374 anti-THRSP antibody WB image

Western blot: 20 µg of SK-BR-3 cell lysate stained with ARG59374 anti-THRSP antibody at 1:2000 dilution.