

## ARG59837 anti-OAF antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes OAF
Tested Reactivity	Ms, Rat
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	OAF
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 28-150 of Human OAF (NP_848602.1).
Conjugation	Un-conjugated
Alternate Names	NS5ATP13TP2; HCV NS5A-transactivated protein 13 target protein 2; Out at first protein homolog

### Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Rat brain and Mouse brain	
Observed Size	30 kDa	

### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.


## Bioinformation

---

Gene Symbol	OAF
Gene Full Name	out at first homolog
Calculated Mw	31 kDa

## Images

---



ARG59837 anti-OAF antibody WB image

Western blot: 25 µg of Rat brain and Mouse brain lysates stained with ARG59837 anti-OAF antibody at 1:3000 dilution.