

ARG62364 anti-Astrocytes antibody [ASTRO1 (10E4/R5)]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [ASTRO1 (10E4/R5)] recognizes Astrocytes
Tested Reactivity	Hu
Tested Application	IHC-P
Host	Mouse
Clonality	Monoclonal
Clone	ASTRO1 (10E4/R5)
Isotype	IgM
Target Name	Astrocytes
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	IHC-P	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Cerebellum	

Properties

Form	Liquid
Purification	Purified Antibody
Buffer	1X PBS and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Astrocyte antibodies ; Astrocyte ELISA Kits ; Astrocyte Duos / Panels ; Astrocyte recombinant proteins ; Anti-Mouse IgM secondary antibodies ;
-----------	---

Related news:

[Astrocyte-to-neuron conversion for Parkinson's disease treatment](#)

Research Area

Cell Biology and Cellular Response antibody; Developmental Biology antibody; Gene Regulation antibody