

**ARG62372**  
anti-Blood Group antigen B antibody [B460]Package: 100 µl  
Store at: -20°C

### Summary

---

Product Description	Mouse Monoclonal antibody [B460] recognizes Blood Group antigen B
Tested Reactivity	Hu
Tested Application	IHC-Fr, IHC-P
Specificity	stains erythrocytes and vascular epithelium of blood group B controls.
Host	Mouse
Clonality	Monoclonal
Clone	B460
Isotype	IgM
Target Name	Blood Group antigen B
Species	Human
Immunogen	Blood Group B erythrocytes
Conjugation	Un-conjugated

### Application Instructions

---

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
Positive Control	Human Tissue (Blood Group B or AB).

### Properties

Form	Liquid
Buffer	