

ARG62390 anti-beta Catenin antibody

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes beta Catenin
Tested Reactivity	Hu, Ms, Dog
Tested Application	ICC/IF, IHC-P, IP, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	beta Catenin
Species	Human
Immunogen	Human beta Catenin.
Conjugation	Un-conjugated
Alternate Names	EVR7; CTNNB; MRD19; armadillo

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A purified
Buffer	