

ARG62440
anti-CDw60 antibody [Jones]

Package: 100 µl

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [Jones] recognizes CDw60
Tested Reactivity	Rat
Tested Application	FACS, ICC/IF
Host	Mouse
Clonality	Monoclonal
Clone	Jones
Isotype	IgM
Target Name	CDw60
Species	Rat
Immunogen	Rat renal epithelia
Conjugation	Un-conjugated

Application Instructions

Application Note	IF: 1: 50-1: 500 FACS: 1µg for 106 cells * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Buffer	PBS, < 0.1% Sodium azide and 0.1% Gelatin
Preservative	< 0.1% Sodium azide
Stabilizer	0.1% Gelatin
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System antibody
---------------	------------------------