

ARG62480
anti-ERCC1 antibody [3H11]Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [3H11] recognizes ERCC1
Tested Reactivity	Hu
Tested Application	IP, WB
Host	Mouse
Clonality	Monoclonal
Clone	3H11
Isotype	IgG2a
Target Name	ERCC1
Species	Human
Immunogen	His-tagged recombinant human ERCC1 protein.
Conjugation	Un-conjugated
Alternate Names	DNA excision repair protein ERCC-1; RAD10; COFS4; UV20

Application Instructions

Application Note	WB: 1/375 - 1/750 IP: 2 µg/mg of lysate * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A purified
Buffer	