

ARG62525 anti-ITK antibody [2F12]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [2F12] recognizes ITK
Tested Reactivity	Hu
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Clone	2F12
Isotype	IgG1
Target Name	ITK
Immunogen	amino-terminal peptide of Itk
Conjugation	Un-conjugated
Alternate Names	LPFS1; Interleukin-2-inducible T-cell kinase; T-cell-specific kinase; Tyrosine-protein kinase ITK/TSK; PSCTK2; EMT; Kinase EMT; Tyrosine-protein kinase Lyk; LYK; IL-2-inducible T-cell kinase; EC 2.7.10.2

Application Instructions

Application Note	WB: 1/1000 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Buffer	