

**ARG62536**  
anti-Firefly Luciferase antibody [LCF01]Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [LCF01] recognizes Firefly Luciferase
Tested Reactivity	Firef
Tested Application	CCle, WB
Host	Mouse
Clonality	Monoclonal
Clone	LCF01
Isotype	IgG1
Target Name	Firefly Luciferase
Species	Firefly
Immunogen	Luciferase protein from Photinus pyralis
Conjugation	Un-conjugated

### Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

### Properties

Form	Liquid
Purification	Protein G affinity purified
Buffer	10mM PBS (pH 7.4), 0.2% BSA and 0.09% Sodium azide
Preservative	0.09% Sodium azide
Stabilizer	0.2% BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Symbol	Phum_PHUM561850
Gene Full Name	firefly luciferase, putative

