

ARG62537 anti-Macrophage antibody [HAM56]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [HAM56] recognizes Macrophage
Tested Reactivity	Hu
Tested Application	ICC/IF, IHC
Host	Mouse
Clonality	Monoclonal
Clone	HAM56
Isotype	IgM, kappa
Target Name	Macrophage
Species	Human
Immunogen	Sonicated human alveolar macrophage.
Conjugation	Un-conjugated

Application Instructions

Application Note	ICC/IF: 1/20 - 1/50 IHC: 1/20 - 1/50 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

Properties

Form	Liquid
Buffer	