

ARG62542 anti-Medroxyprogesterone Acetate antibody [MP01(1F5)]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rat Monoclonal antibody [MP01 (1F5)] recognizes Medroxyprogesterone Acetate
Tested Reactivity	Other
Tested Application	ELISA, RIA
Host	Rat
Clonality	Monoclonal
Clone	MP01 (1F5)
Isotype	IgG1
Target Name	Medroxyprogesterone Acetate
Species	Others
Immunogen	Medroxyprogesterone Acetate (MPA) -3-carboxymethyl oxime BSA.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	RIA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified Antibody
Buffer	1X PBS and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

