

ARG62588
anti-PTHrP antibody [3D1.1]Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [3D1.1] recognizes PTHR1
Tested Reactivity	Hu
Tested Application	IHC, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	3D1.1
Isotype	IgG1, kappa
Target Name	PTHrP
Species	Human
Immunogen	Synthetic peptide corresponding to amino acids 146-169 of the human PTH/PTHrP receptor (RRCDRNGSWELVPGHNRTWANYSE).
Conjugation	Un-conjugated
Alternate Names	Parathyroid hormone/parathyroid hormone-related peptide receptor; PTH/PTHrP type I receptor; Parathyroid hormone 1 receptor; PTHR1; PTH/PTHr receptor; PFE; PTHR; PTH1 receptor

Application Instructions

Application Note	WB: 0.5-2 µg/ml IHC: 1/10 - 1/500 IHC-P: 5 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein G purified
Buffer	