

ARG62600 anti-Plasma Cell antibody [LIV3G11]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [LIV3G11] recognizes Plasma Cell
Tested Reactivity	Hu
Species Does Not React With	Rat
Tested Application	IHC-P
Host	Mouse
Clonality	Monoclonal
Clone	LIV3G11
Isotype	IgG2a
Target Name	Plasma Cell
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	IHC-P	1:400
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Tonsil or Lymph node	

Properties

Form	Liquid
Purification	Purified Antibody
Buffer	1X PBS and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody; Immune System antibody
---------------	---