

ARG62609 anti-PKC antibody [MC5]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [MC5] recognizes PKC
Tested Reactivity	Hu, Ms, Rat
Tested Application	ELISA, FACS, ICC/IF, IHC-Fr, IHC-P, IP, RIA, WB
Host	Mouse
Clonality	Monoclonal
Clone	MC5
Isotype	IgG2a
Target Name	PKC
Species	Bovine
Immunogen	Purified bovine brain protein kinase C
Conjugation	Un-conjugated
Alternate Names	PKC-alpha; EC 2.7.11.13; AAG6; PKCA; PRKACA; PKC-A; Protein kinase C alpha type

Application Instructions

Application Note	FACS: 1µg for 106 cells IHC-Fr: 1/500 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A purified
Buffer	