

## ARG62610 anti-PNRA antibody [B461]

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [B461] recognizes PNRA
Tested Reactivity	Hu, Rat
Tested Application	IHC-Fr, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	B461
Isotype	IgG2b
Target Name	PNRA
Immunogen	Renal cortical tissue extract
Conjugation	Un-conjugated

### Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--------------------------------------------------------------------------------------------------------------------------------------------

### Properties

Form	Liquid
Purification	Protein A affinity purified
Buffer	PBS (pH 7.4), 0.2% BSA and 0.09% Sodium azide
Preservative	0.09% Sodium azide
Stabilizer	0.2% BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Cancer antibody
---------------	-----------------