

ARG62618 anti-Retinoic Acid Receptor alpha antibody [763]

Package: 100 µl
Store at: -20°C

Summary


Product Description	Mouse Monoclonal antibody [763] recognizes Retinoic Acid Receptor alpha
Tested Reactivity	Hu, Ms, Rat
Predict Reactivity	Bov, Dog, Gpig
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Clone	763
Isotype	IgG1
Target Name	Retinoic Acid Receptor alpha
Species	Human
Immunogen	Synthetic peptide corresponding to the N-terminal region of RARA conjugated to KLH
Conjugation	Un-conjugated
Alternate Names	Nuclear receptor subfamily 1 group B member 1; RAR; RAR-alpha; Retinoic acid receptor alpha; NR1B1

Application Instructions

Application Note	WB: 1/1000 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein G purified
Buffer	



ARG62618 anti-Retinoic Acid Receptor alpha antibody [763] WB image

Western blot: Rat hippocampal lysate showing specific immunolabeling of the ~48 kDa RAR-alpha protein stained with ARG62618 anti-Retinoic Acid Receptor alpha antibody [763].