

## ARG62689 anti-Blood Group A1B antigen antibody [HE-24]

Package: 1 ml  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [HE-24] recognizes Blood Group A1B antigen
Tested Reactivity	Hu
Tested Application	Agg
Specificity	The clone HE-24 distinguishes blood group A <sub>1</sub> B from A <sub>2</sub> B. The specificity of the antibody HE-24 was verified on >1000 samples of blood.
Host	Mouse
Clonality	Monoclonal
Clone	HE-24
Isotype	IgM
Target Name	Blood Group A1B antigen
Immunogen	Glycoprotein fraction isolated from saliva of secretors with blood group A.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	Agg	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification Note	The hybridoma culture supernatant is 4 x concentrated by ultrafiltration using 100 kDa-cut off membrane.
Buffer	Hybridoma culture supernatant
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Cell Biology and Cellular Response antibody
---------------	---

Clone	Agglutination of tested blood groups				
	A <sub>1</sub>	B	A <sub>1</sub> B	A <sub>2</sub> B	0
HE-193 (anti-A)	+	-	+	+	-
HE-195 (anti-A)	+	-	+	+	-
HEB-29 (anti-B)	-	+	+	+	-
HE-24 (anti-A <sub>1</sub> B)	+	-	+	-	-

ARG62689 anti-Blood Group A1B antigen antibody [HE-24]  
Agglutination image

Agglutination of particular blood groups using ARG62689 anti-Blood Group A1B antigen antibody [HE-24] (anti-blood group A1B).

---