

ARG66220 anti-Aflatoxin B1 + Aflatoxin B2 antibody [BCAB1311]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [BCAB1311] recognizes Aflatoxin B1 + Aflatoxin B2
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Cross-Reactivity: Aflatoxin B2 (100%), Aflatoxin B1 (65%), Aflatoxin G2 (11%), Aflatoxin G1 (6%)
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1311
Isotype	IgG1
Target Name	Aflatoxin B1 + Aflatoxin B2
Species	Others
Immunogen	HAS-conjugated Aflatoxin B2.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95%. Purity is tested by electrophoresis.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide.
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Calculated Mw

150 kDa