

**ARG66221**  
anti-Zearalenone antibody [BCAB1330]Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [BCAB1330] recognizes Zearalenone
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Cross-Reactivity: Zearalenone (100%), Zearalanone (60%), alpha-Zearalanol (40%), alpha-Zearalenol (11%), beta-Zearalanol (3%), beta-Zearalenol (1%)
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1330
Isotype	IgG1
Target Name	Zearalenone
Species	Others
Immunogen	BSA-conjugated Zearalenone
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95%. Purity is tested by electrophoresis.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide.
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.