

## ARG66303 Rabbit anti-Duck IgG (H+L) antibody (HRP)

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Rabbit Polyclonal antibody recognizes Duck IgG (H+L)
Tested Reactivity	Duck
Tested Application	Dot, ELISA, ICC, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Duck
Immunogen	Duck IgG
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	Dot	Assay-dependent
	ELISA	Assay-dependent
	ICC	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Purified.
Buffer	PBS (pH 7.2)
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.