

**ARG66743**  
anti-TAP tag antibodyPackage: 100 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody recognizes TAP tag
Tested Reactivity	Other
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Target Name	TAP tag
Species	Others
Immunogen	Recombinant Protein of TAP tag.
Epitope	Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	WB	1:500 - 1:10000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4), 0.02% Sodium azide, 50% Glycerol and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol and 0.5% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



TAP recombinant protein

ARG66743 anti-TAP tag antibody WB image

Western blot: TAP recombinant protein stained with ARG66743 anti-TAP tag antibody at 1:5000 dilution.