

ARG66744
anti-TAP tag antibody (HRP)Package: 100 µg
Store at: 4°C

Summary

Product Description	HRP-conjugated Mouse Monoclonal antibody recognizes TAP tag
Tested Reactivity	Other
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Target Name	TAP tag
Species	Others
Immunogen	Recombinant Protein of TAP tag.
Epitope	Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag.
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	WB	1:5000 - 1:10000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.