

## ARG80372 Human C1q antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

|                     |   |
|---------------------|---|
| Product Description | ARG80372 Human C1q antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human C1q Ab in serum and plasma (EDTA, heparin, citrate). |
| Tested Reactivity   | Hu  |
| Tested Application  | ELISA   |
| Target Name         | C1q antibody  |
| Sensitivity         | 0.5 U/ml  |
| Sample Type         | Serum and plasma (EDTA, heparin, citrate).  |
| Standard Range      | 6.3 - 100 U/ml  |
| Sample Volume       | 10 µl   |

### Application Instructions

|            |                        |
|------------|------------------------|
| Assay Time | 30, 15, 15, 5 min (RT) |
|------------|------------------------|

### Properties

|                     |  |
|---------------------|--|
| Form                | 96 well  |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.   |

### Bioinformation

|               |   |
|---------------|---|
| Highlight     | Related products:<br><a href="#">Complement component antibodies:</a> <a href="#">Complement component ELISA Kits:</a> <a href="#">Complement component Duos / Panels:</a><br>New ELISA data calculation tool:<br><a href="#">Simplify the ELISA analysis by GainData</a> |
| Research Area | Immune System kit; Neuroscience kit   |