

ARG80534 Human Chlamydia trachomatis IgG antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80534 Human Chlamydia trachomatis IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of IgG class antibodies against Chlamydia trachomatis in serum and plasma (heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Chlamydia trachomatis IgG antibody
Sensitivity	diagn. >95%
Sample Type	Serum and plasma (heparin, citrate).
Standard Range	Cut-off
Sample Volume	10 µl

Application Instructions

Assay Time 1 h (37°C), 30 min (RT), 15 min (RT/dark)

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area Immune System kit; Microbiology and Infectious Disease kit