

ARG80821 Human Allergen-related IgG4 antibody screen ELISA Kit (20 Food Antigenen)

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG80821-001	Microtiter plate	12 strips X 8 wells (for 4 patients).	4°C
ARG80821-002	Standards (0.35, 0.70, 3.5, 17.5 U/ml IgG4 antibodies against f01)	4 X 0.5 ml (ready to use)	4°C
ARG80821-003	Sample dilution buffer	40 ml (ready to use)	4°C
ARG80821-004	Mouse anti human IgG4-AP-conjugate	15 ml (ready to use)	4°C
ARG80821-005	10X Wash buffer	60 ml	4°C
ARG80821-006	PNPP substrate	15 ml (ready to use)	4°C (Protect from light)
ARG80821-007	STOP solution	15 ml (ready to use)	4°C

Summary

Product Description	ARG80821 Human Allergen-related IgG4 antibody screen ELISA Kit (20 Food Antigenen) (20 Food Antigenen) is an Enzyme Immunoassay kit for screening Human's Allergen-related IgG4 antibody for 20 Food Antigenen in serum and plasma (EDTA, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Allergen-related IgG4 antibody screen
Conjugation	HRP
Conjugation Note	PNPP at 405 nm
Sample Type	Serum and plasma (EDTA, heparin).
Standard Range	0.35 - 17.5 U/ml
Sample Volume	25 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note

For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area

Immune System kit