

## ARG80991 Dog CPV / Canine Parvo Virus IgM antibody ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	ARG80991 Dog CPV / Canine Parvo Virus IgM antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Dog's CPV / Canine Parvo Virus IgM antibody in serum and plasma.
Tested Reactivity	Dog
Tested Application	ELISA
Target Name	Dog CPV / Canine Parvo Virus IgM
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum and plasma.
Sample Volume	100 µl

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---