

ARG80992 Dog CPV / Canine Parvo Virus antigen ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80992 Dog CPV / Canine Parvo Virus antigen ELISA Kit is an Enzyme Immunoassay kit for the determination of Dog CPV / Canine Parvo Virus antigen in faeces (supernatant of dilution).
Tested Reactivity	Virus
Tested Application	ELISA
Target Name	Dog CPV / Canine Parvo Virus antigen
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Faeces (supernatant of dilution).
Sample Volume	50 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---