

## ARG80998 Dog CRV / Canine Rabies Virus IgG antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80998 Dog CRV / Canine Rabies Virus IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Canine Rabies Virus IgG antibody in dog's serum and plasma.
Tested Reactivity	Dog
Tested Application	ELISA
Target Name	Dog CRV / Canine Rabies Virus IgG
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum and plasma.
Sample Volume	100 µl

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---