

## ARG81132 Nitrotyrosine ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG81132 Nitrotyrosine ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Nitrotyrosine in plasma, serum, cell lysates, urine, and other sample matrices.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Nitrotyrosine
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm
Sensitivity	50 nM
Sample Type	Plasma, serum, cell lysates, urine, and other sample matrices.
Standard Range	62.5 - 8000 nM
Sample Volume	50 µl

### Application Instructions


Assay Time	~ 2 hours
------------	-----------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">Nitrotyrosine antibodies</a> ; <a href="#">Nitrotyrosine ELISA Kits</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
-----------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



ARG81132 Nitrotyrosine ELISA Kit standard curve image

ARG81132 Nitrotyrosine ELISA Kit results of a typical standard run with optical density reading at 450 nm.