

ARG81202 Cholic Acid ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG81202 Cholic Acid ELISA Kit is an Enzyme Immunoassay kit for the quantification of Cholic Acid in plasma, serum, urine, feces, cell and tissue lysates.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Cholic Acid
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	0.4 µM
Sample Type	Plasma, serum, urine, feces, and cell and tissue lysates.
Standard Range	0.39 - 25 µM
Sample Volume	50 µl

Application Instructions


Assay Time	~ 3 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG81202 Cholic Acid ELISA Kit standard curve image

ARG81202 Cholic Acid ELISA Kit results of a typical standard run with optical density reading at 450 nm.