

ARG81221
CML-HDL ELISA KitPackage: 96 wells
Store at: 4°C, -80°C

Summary

Product Description	ARG81221 CML-HDL ELISA Kit is an Enzyme Immunoassay kit for the quantification of CML-HDL in plasma or serum.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	CML-HDL
Conjugation Note	Read at 450 nm.
Sensitivity	< 1 ng/ml
Sample Type	Plasma or serum
Standard Range	0.31 - 20 ng/ml
Sample Volume	100 µl

Application Instructions


Assay Time	4.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -80°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG81221 CML-HDL ELISA Kit standard curve image

ARG81221 CML-HDL ELISA Kit results of a typical standard run with optical density reading at 450 nm.