

ARG81344 HBV surface antigen / HBsAg ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG81344 Hepatitis B surface antigen / HBs Ag ELISA Kit is an Enzyme Immunoassay kit for the quantification of Hepatitis B surface antigen / HBs Ag in purified virus or unpurified viral cell culture supernatant.
Tested Reactivity	HBV
Tested Application	ELISA
Target Name	HBV surface antigen / HBsAg
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	~ 1 ng/ml
Sample Type	Purified virus or unpurified viral cell culture supernatant.
Standard Range	1.56 - 100 ng/ml
Sample Volume	100 µl

Application Instructions


Assay Time	~ 4.5 hours
------------	-------------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: HBV antibodies ; HBV ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	---



ARG81344 HBV surface antigen / HBsAg ELISA Kit standard curve image

ARG81344 HBV surface antigen / HBsAg ELISA Kit results of a typical standard run with optical density reading at 450 nm.