


ARG81358 Mouse/Rat CGRP ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81358 Mouse/Rat CGRP ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse/Rat CGRP in plasma.
Tested Reactivity	Ms, Rat
Tested Application	ELISA
Specificity	

Images



ARG81358 Mouse/Rat CGRP ELISA Kit standard curve image

ARG81358 Mouse/Rat CGRP ELISA Kit results of a typical standard run with optical density reading at 450 nm.