

ARG81362 HCV Core antigen / HCcAg ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG81362 HCV core antigen ELISA Kit is an Enzyme Immunoassay kit for the quantification of HCV core antigen in purified virus or unpurified viral supernatant.
Tested Reactivity	HCV
Tested Application	ELISA
Target Name	HCV Core antigen / HCcAg
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	1 ng/ml
Sample Type	Purified virus or unpurified viral supernatant.
Standard Range	1.56 - 100 ng/ml
Sample Volume	100 µl

Application Instructions


Assay Time	~ 4.5 hours
------------	-------------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: HCV antibodies ; HCV ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	---



ARG81362 HCV Core antigen / HCcAg ELISA Kit standard curve image

ARG81362 HCV Core antigen / HCcAg ELISA Kit results of a typical standard run with optical density reading at 450 nm.