

## ARG81628 Triglyceride Assay Kit

Package: 200 tests  
Store at: -20°C

### Summary


Product Description	ARG81628 Triglyceride Assay Kit is a detection kit for the quantification of Triglyceride in serum and plasma.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Triglyceride
Conjugation Note	Read at 570 nm.
Sensitivity	0.01 mmol/L (0.88 mg/dL)
Sample Type	Serum and plasma.
Standard Range	0.3 - 1 mmol/L
Sample Volume	10 µl
Alternate Names	TG; triacylglycerol; TAG; or triacylglyceride

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	~ 1 hour

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG81628 Triglyceride Assay Kit standard curve image

ARG81628 Triglyceride Assay Kit results of a typical standard run with optical density reading at 570 nm.