

ARG81629 HDL and LDL / VLDL Assay Kit

Package: 100 tests
Store at: -20°C

Summary

| | |
|---------------------|---|
| Product Description | ARG81629 HDL and LDL / VLDL Assay Kit is a detection kit for the quantification of HDL and LDL / VLDL in serum. |
| Tested Reactivity | Other |
| Tested Application | FuncSt |
| Target Name | HDL and LDL / VLDL |
| Conjugation Note | Read at 340 nm. |
| Sensitivity | 5 mg/dL |
| Sample Type | Serum. |
| Sample Volume | 20 µl |

Application Instructions

| | |
|------------------|--|
| Application Note | Please note that this kit does not include a microplate. |
| Assay Time | ~ 1 hour |

Properties

| | |
|---------------------|---|
| Form | Liquid |
| Storage instruction | Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |