

ARG82020 Tissue Glucose Assay Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary


Product Description	ARG82020 Tissue Glucose Assay Kit is a detection kit for the quantification of Glucose Content in tissue extracts and cell lysate.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Tissue Glucose
Conjugation Note	Read at 505 nm.
Sensitivity	0.1 mmol/l
Detection Range	0.1 - 10 mmol/l
Sample Type	Tissue extracts and cell lysate.
Sample Volume	20 µl

Application Instructions

Assay Time	15 min
------------	--------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82020 Tissue Glucose Assay Kit typical data demonstration image

ARG82020 Tissue Glucose Assay Kit results of a typical data with optical density reading at 505 nm.