

ARG82055 Human eNOS ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	Human eNOS ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human eNOS in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	eNOS
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	39 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	78 - 5000 pg/ml
Sample Volume	100 µl
Alternate Names	Constitutive NOS; NOS type III; Nitric oxide synthase, endothelial; Endothelial NOS; eNOS; EC-NOS; NOSIII; cNOS; EC 1.14.13.39; ECNOS

Application Instructions

Assay Time	4.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	NOS3
Gene Full Name	nitric oxide synthase 3 (endothelial cell)
Background	Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. Nitric oxide is synthesized from L-arginine by nitric oxide synthases. Variations in this gene are associated with susceptibility to coronary spasm. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009]
Function	Produces nitric oxide (NO) which is implicated in vascular smooth muscle relaxation through a cGMP-mediated signal transduction pathway. NO mediates vascular endothelial growth factor (VEGF)-induced angiogenesis in coronary vessels and promotes blood clotting through the activation of platelets.

Isoform eNOS13C: Lacks eNOS activity, dominant-negative form that may down-regulate eNOS activity by forming heterodimers with isoform 1. [UniProt]

Highlight

Related products:
[eNOS antibodies](#); [eNOS ELISA Kits](#); [eNOS Duos / Panels](#);
New ELISA data calculation tool:
[Simplify the ELISA analysis by GainData](#)


PTM

Phosphorylation by AMPK at Ser-1177 in the presence of Ca(2+)-calmodulin (CaM) activates activity. In absence of Ca(2+)-calmodulin, AMPK also phosphorylates Thr-495, resulting in inhibition of activity (By similarity). Phosphorylation of Ser-114 by CDK5 reduces activity. [UniProt]

Cellular Localization

Cell membrane. Membrane, caveola. Cytoplasm, cytoskeleton. Golgi apparatus. Note=Specifically associates with actin cytoskeleton in the G2 phase of the cell cycle; which is favored by interaction with NOSIP and results in a reduced enzymatic activity. [UniProt]

Images



ARG82055 Human eNOS ELISA Kit standard curve image

ARG82055 Human eNOS ELISA Kit results of a typical standard run with optical density reading at 450 nm.
