

ARG82096 Mouse Uteroglobin ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG82096-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82096-002	Standard	2 X 10 ng/vial	4°C
ARG82096-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82096-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82096-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82096-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82096-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82096-008	25X Wash buffer	20 ml	4°C
ARG82096-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82096-010	STOP solution	10 ml (Ready to use)	4°C
ARG82096-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG82096 Mouse Uteroglobin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse Uteroglobin in serum, plasma (heparin, EDTA), urine and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	Uteroglobin
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	2.35 pg/ml
Sample Type	Serum, plasma (heparin, EDTA), urine and cell culture supernatants.
Standard Range	4.7 - 300 pg/ml
Sample Volume	100 µl

Precision	Intra-Assay CV: 5.2%; Inter-Assay CV: 5.9%
Alternate Names	Secretoglobin family 1A member 1; Urine protein 1; Urinary protein 1; Clara cell phospholipid-binding protein; UGB; CCPBP; CCSP; UP-1; CC10; Clara cells 10 kDa secretory protein; UP1; Uteroglobin; CC16

Application Instructions

Assay Time	~ 5 hours
------------	-----------


Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	SCGB1A1
Gene Full Name	secretoglobin, family 1A, member 1 (uteroglobin)
Background	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [provided by RefSeq, May 2010]
Function	Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2. [UniProt]
Highlight	Related products: Uteroglobin antibodies ; Uteroglobin ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Cellular Localization	Secreted. [UniProt]

Images



ARG82096 Mouse Uteroglobin ELISA Kit standard curve image

ARG82096 Mouse Uteroglobin ELISA Kit results of a typical standard run with optical density reading at 450 nm.