

ARG82141 Bilirubin Assay Kit (Colorimetric)

Package: 180 tests
Store at: 4°C

Summary


Product Description	ARG82141 Bilirubin Assay Kit (Colorimetric) is a detection kit for the quantification of Bilirubin in serum.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Bilirubin
Conjugation Note	Read at 530 nm
Sensitivity	0.16 mg/dL
Detection Range	1 - 10 mg/dL
Sample Type	serum
Sample Volume	50 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	10 min

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82141 Bilirubin Assay Kit (Colorimetric) typical data demonstration image

ARG82141 Bilirubin Assay Kit (Colorimetric) results of a typical data with optical density reading at 530 nm.