

ARG82146 Coenzyme A Assay Kit (Colorimetric)

Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82146 Coenzyme A Assay Kit (Colorimetric) is a detection kit for the quantification of Coenzyme A in Biological samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Coenzyme A
Conjugation Note	Read at 570 nm (colorimetric), or 530/585 nm (fluorometric)
Sensitivity	3 μ M
Detection Range	5 - 1000 μ M (colorimetric); 3 - 100 μ M (fluorometric)
Sample Type	Biological samples
Sample Volume	10 μ l

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	60 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82146 Coenzyme A Assay Kit (Colorimetric) typical data demonstration image

ARG82146 Coenzyme A Assay Kit (Colorimetric) results of a typical data with optical density reading at 570 nm.