

ARG82156 Fructose Assay Kit (Colorimetric)

Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82156 Fructose Assay Kit (Colorimetric) is a detection kit for the quantification of Fructose in serum, plasma, urine, saliva, milk, culture medium, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Fructose
Conjugation Note	Read at 565 nm
Sensitivity	12 µM
Detection Range	12 - 1000 µM
Sample Type	serum, plasma, urine, saliva, milk, culture medium, food and beverage
Sample Volume	20 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	70 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82156 Fructose Assay Kit (Colorimetric) typical data demonstration image

ARG82156 Fructose Assay Kit (Colorimetric) results of a typical data with optical density reading at 565 nm.