

## ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric)

Package: 100 tests  
Store at: -20°C

### Summary


|                     |  |
|---------------------|--|
| Product Description | ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) is a detection kit for the quantification of Glucose 6 phosphate in serum, plasma, tissue and cell culture supernatants. |
| Tested Reactivity   | Hu, Ms, Rat, All   |
| Tested Application  | FuncSt   |
| Target Name         | Glucose 6 phosphate  |
| Conjugation Note    | Read at 460 nm   |
| Sensitivity         | 10 µM  |
| Detection Range     | 10 - 1000 µM   |
| Sample Type         | serum, plasma, tissue and cell culture supernatants  |
| Sample Volume       | 20 µl  |

### Application Instructions

|                  |  |
|------------------|--|
| Application Note | Please note that this kit does not include a microplate. |
| Assay Time       | 20 min   |

### Properties

|                     |   |
|---------------------|---|
| Form                | Liquid  |
| Storage instruction | Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.  |



ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) typical data demonstration image

ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) results of a typical data with optical density reading at 460 nm.