

## ARG82170 Indole Assay Kit (Colorimetric)

Package: 100 tests  
Store at: 4°C

### Summary


Product Description	ARG82170 Indole Assay Kit (Colorimetric) is a detection kit for the quantification of Indole in bacterial growth medium.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Indole
Conjugation Note	Read at 565 nm
Sensitivity	3 µM
Detection Range	3 - 100 µM
Sample Type	bacterial growth medium
Sample Volume	100 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	< 5 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82170 Indole Assay Kit (Colorimetric) typical data demonstration image

ARG82170 Indole Assay Kit (Colorimetric) results of a typical data with optical density reading at 565 nm.