

ARG82186 Oxalate Assay Kit (Colorimetric)

Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82186 Oxalate Assay Kit (Colorimetric) is a detection kit for the quantification of Oxalate in urine, animal and plant tissue samples..
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Oxalate
Conjugation Note	Read at 595 nm
Sensitivity	20 µM
Detection Range	20 - 1500 µM
Sample Type	urine, animal and plant tissue samples.
Sample Volume	10 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	10 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82186 Oxalate Assay Kit (Colorimetric) typical data demonstration image

ARG82186 Oxalate Assay Kit (Colorimetric) results of a typical data with optical density reading at 595 nm.