

## ARG82190 Phenylalanine Assay Kit (Fluorometric)

Package: 100 tests  
Store at: -20°C

### Summary


Product Description	ARG82190 Phenylalanine Assay Kit (Fluorometric) is a detection kit for the quantification of Phenylalanine in serum and urine.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Phenylalanine
Conjugation Note	Read at 530/595 nm (fluorometric)
Sensitivity	2 µM
Detection Range	2 - 300 µM
Sample Type	serum and urine
Sample Volume	10 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82190 Phenylalanine Assay Kit (Fluorometric) typical data demonstration image

ARG82190 Phenylalanine Assay Kit (Fluorometric) results of a typical data.