

## ARG82195 Sulfate Assay Kit

Package: 200 tests  
Store at: 4°C

### Summary


Product Description	ARG82195 Sulfate Assay Kit is a detection kit for the quantification of Sulfate in serum, urine and food.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Sulfate
Conjugation Note	Read at 600 nm
Sensitivity	20 µM
Detection Range	0.02 mM (0.19 mg/dL) - 2 mM (19.2 mg/dL)
Sample Type	serum, urine and food
Sample Volume	200 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	5 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82195 Sulfate Assay Kit typical data demonstration image

ARG82195 Sulfate Assay Kit results of a typical data with optical density reading at 600 nm.